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Experimental Report on H&E Staining

1. Instruments and key reagents

1.1 Instruments

Instrument	Manufacture	Specifications/Model
Tissue processor	Changzhou Zhongwei Electronics Co., Ltd	TSJ-SD
Tissue embedder	Changzhou Zhongwei Electronics Co., Ltd	BMJ-A
Microtome	ThermoFisher Scientific	SHANDON FINESSE 325
Freezing table	Wuhan Junjie Electronics Co., Ltd	JB-L5
Water bath - Slide drier	Changzhou Zhongwei Electronics Co., Ltd	PHY-III
Slide	Hubei BIOSSCI Biotech Co., Ltd	BP0510
Upright microscope	Olympus	CX-31
Digital scanner	HAMAMATSU PHOTONICS	NanoZoomer®S360

1.2 Key reagents

Reagent	Manufacture	Specifications/Model
Ethanol	Sinopharm	100092683
Clearer	Wuhan Tongsheng Technology Development Co., Ltd	
Neutral balsam	Wuhan Tongsheng Technology Development Co., Ltd	
H&E staining kit	Hubei BIOSSCI Biotech Co., Ltd	BP0211

2. Procedures

2.1 Deparaffinization and rehydration.

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2.1.1 Tissue sections were immersed in clearer for 10min. Repeat this step two times, gently shaking off excess liquid between each step.

2.1.2 Tissue sections were immersed in progressively more dilute ethanol solutions and ultimately immersed in distilled water to rehydrate the tissue: Absolute ethanol for 5min, 95% ethanol for 5min, 85% ethanol for 5min, 75% ethanol for 5min. Rinsing with distilled water for 1min.

2.2 Sections were stained with hematoxylin solution (Harris) for 4 minutes and then washed with tap water for 2 minutes until no excess dye comes out of the sections.

2.3 Sections were differentiated with 0.8% hydrochloric acid alcohol for 2 seconds and washed with tap water. Sections can also be treated with lithium carbonate solution to be bluer and washed with tap water for 2 minutes.

2.4 Sections were stained with eosin solution (alcohol soluble) for 20 seconds without water washing. Then, sections were treated with 95% ethanol for 5 minutes and dehydrated with absolute ethanol I and II for 2 minutes.

2.5 The tissue sections were transparent with xylene and then mounted with neutral balsam. Examination with microscope.

3. Results

The nucleus is blue and the cytoplasm is pink to red.

4. Attentions

4.1 It is important that paraffin sections need to be fully dewaxed.

4.2 If nuclei are too dark, it can be differentiated again.

5. Technical support

The technical support of H&E staining experiment is provided by Hubei BIOSSCI Biotech Co., Ltd (Wuhan Changyan Pathology technology Co., Ltd).

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