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Experimental Report on Picro Sirius Red Staining

1. Instruments and key reagents

1.1 Instruments

Instrument	Manufacture	Specifications/Model
Tissue processor	Changzhou Zhongwei Electronics Co., Ltd	TSJ-SD
Tissue embedder	Changzhou Zhongwei Electronics Co., Ltd	BMJ-A
Microtome	ThermoFisher Scientific	SHANDON FINESSE 325
Freezing table	Wuhan Junjie Electronics Co., Ltd	JB-L5
Water bath - Slide drier	Changzhou Zhongwei Electronics Co., Ltd	PHY-III
Slide	Hubei BIOSSCI Biotech Co., Ltd	BP0510
Upright microscope	Olympus	CX-31
Digital scanner	HAMAMATSU PHOTONICS	NanoZoomer®S360

1.2 Key reagents

Reagent	Manufacture	Specifications/Model
Ethanol	Sinopharm	100092683
Clearer	Wuhan Tongsheng Technology Development Co., Ltd	
Neutral balsam	Wuhan Tongsheng Technology Development Co., Ltd	
Picro Sirius Red staining kit	Hubei BIOSSCI Biotech Co., Ltd	

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2.1 Tissues were fixed with neutral formaldehyde fix solution.

2.2 Deparaffinization and rehydration.

2.2.1 Tissue sections were immersed in clearer for 10min. Repeat this step two times, gently shaking off excess liquid between each step.

2.2.2 Tissue sections were immersed in progressively more dilute ethanol solutions and ultimately immersed in distilled water to rehydrate the tissue: Absolute ethanol for 5min, 95% ethanol for 5min, 85% ethanol for 5min, 75% ethanol for 5min. Rinsing with distilled water for 1min.

2.3 The tissue sections were stained with celestine blue solution for 5-10 minutes and washed with distilled water. (Harris hematoxylin can also be used instead)

2.4 The tissue sections were stained with saturated Picro Sirius Red solution for 15-30 minutes.

2.5 The tissue sections were differentiated and dehydrated with absolute ethanol, which required about four cylinders of ethanol for dehydration.

2.6 The tissue sections were transparent with xylene and mounted with neutral balsam.

3. Results

Collagen fibers are red. Nuclei are green and the others are yellow.

4. Attentions

4.1 It is crucial to observe and take photos with a polarization microscope in time to keep bright colors. It can also be observed under an upright microscope, but only red and yellow can be seen.

4.2 There are four types of collagen fibers can be observed under a polarization microscope:

Collagen I: Yellow or red fibers. Collagen fibers are tightly arranged, showing strong birefringence.

Collagen II: It shows weak birefringence and a loose network distribution with

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various colors.

Collagen III: It shows weak birefringence and is a green fine fiber.

Collagen IV: It shows basal lamina with weak birefringence, which is light yellow.

5. Technical support

The technical support of Picro Sirius Red staining experiment is provided by Hubei BIOSSCI Biotech Co., Ltd (Wuhan Changyan Pathology technology Co., Ltd).

Tel: 400 118 0100

Fax: +86-027-87382710

E-mail: support@biossci.com

Website: www.biossci.com