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Experimental Report on Toluidine Blue (Nissl Bodies)

Staining

1. Instruments and key reagents

1.1 Instruments

Instrument	Manufacture	Specifications/Model
Tissue processor	Changzhou Zhongwei Electronics Co., Ltd	TSJ-SD
Tissue embedder	Changzhou Zhongwei Electronics Co., Ltd	BMJ-A
Microtome	ThermoFisher Scientific	SHANDON FINESSE 325
Freezing table	Wuhan Junjie Electronics Co., Ltd	JB-L5
Water bath - Slide drier	Changzhou Zhongwei Electronics Co., Ltd	PHY-III
Slide	Hubei BIOSSCI Biotech Co., Ltd	BP0510
Upright microscope	Olympus	CX-31
Digital scanner	HAMAMATSU PHOTONICS	NanoZoomer®S360

1.2 Key reagents

Reagent	Manufacture	Specifications/Model
Ethanol	Sinopharm	100092683
Clearer	Wuhan Tongsheng Technology Development Co., Ltd	
Neutral balsam	Wuhan Tongsheng Technology Development Co., Ltd	
Toluidine blue solution	Hubei BIOSSCI Biotech Co., Ltd	BP0360

2. Procedures

2.1 Deparaffinization and rehydration.

2.1.1 Tissue sections were immersed in clearer for 10min. Repeat this step two times, gently shaking off excess liquid between each step.

2.1.2 Tissue sections were immersed in progressively more dilute ethanol solutions and ultimately immersed in distilled water to rehydrate the tissue: Absolute ethanol for 5min, 95% ethanol for 5min, 85% ethanol for 5min, 75% ethanol for 5min. Rinsing with distilled water for 1min.

2.2 Tissue sections were placed in preheated (50°C), 1% toluidine blue aqueous solution and sections were stained in incubator at 56°C for 20 minutes. Then, sections were washed with distilled water.

2.3 Tissue sections were differentiated with 95% ethanol or 0.1% glacial acetic acid and controlled under microscope until Nissl bodies were clearly visible.

2.4 Tissue sections were dehydrated rapidly with absolute ethanol.

2.5 Tissue sections were put into clearer and then mounted with neutral balsam. Examination with microscope.

3. Results

Nuclei are pare blue. Nissl bodies are dark blue. The background is light blue or colorless.

4. Attentions

4.1 Tissues should be fixed immediately because the Nissl body is easy to dissolve in vitro, otherwise, it is difficult to be stained.

4.2 The differentiation with 95% ethanol should be carried out quickly. Observe the tissue section with naked eyes until the section is clear and the background is light blue or colorless.

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The technical support of Toluidine Blue (Nissl Bodies) staining experiment is provided by Hubei BIOSSCI Biotech Co., Ltd (Wuhan Changyan Pathology technology Co., Ltd).

Tel: 400 118 0100

Fax: +86-027-87382710

E-mail: support@biossci.com

Website: www.biossci.com